


Innovative solutions for effective diagnostics



NEW PRODUCT


Strep-A Respi-Strip

Rapid test for the detection of the specific antigen of **Group A *Streptococcus*** directly from throat swab specimens.

- Specific detection of **Group A *Streptococcus*** !
- Result after 5 minutes !
- Easy to use and easy to read !



Strep-A Respi-Strip



● Performances ●

	Culture method (n=756)	
Specificity	94.5%	(92 to 96.3%)
Sensitivity	90.6%	(86.2 to 93.7%)
PPV	89.9%	(85.5 to 93.1%)
NPV	94.9%	(92.4 to 96.6%)
Reliability	93.1%	


This evaluation was performed on throat swabs specimens collected from pediatric and adult populations of patients exhibiting symptoms of pharyngitis.

The conventional culture test was used as reference method.

● Ordering Information ●

Product Name	Description	Storage	Code
Strep-A Respi-Strip	25 tests, dilution buffer, extraction buffer	4°C to 30°C	C-1022

Science Park CREALYS
Rue Jean Sonet 4A
5032 Gembloux
BELGIUM



Ph.: +32 (0)81 719 910
Fx.: +32 (0)81 719 919
sales@corisbio.com
www.corisbio.com